

CA800 Chest sling CA900 Transport sling User instructions



CA800



CA900

CA 800 Chest Sling

For standing only



Ensure that the guiding handles and identification label are on the outside and place the sling down the back of the client. Check the sling is equal and square.



Close the Velcro belt and apply hook Velcro to belt. Close belt so it is comfortable.



The additional belt with buckle can now be closed. This belt can be adjusted by pulling the end of strap. Attach straps to hoist, the shortest loop will achieve a more upright standing position.

The CA800 Chest Sling is suitable for those clients who have a degree of weight bearing ability. It is designed to assist the client to stand.

CA 900 Transport Sling

For seat to seat transfer only



Apply the sling in the same way as CA 800.

Raise and support the clients leg one at a time and carefully feed the leg strap underneath and up between both legs.

Ensure that the sling is not twisted or creased under the thighs and check that both leg straps are equal in length. Cross the leg straps

Apply the leg straps on longest loop and the shoulder on a suitable loop to determine the clients posture. The shorter the shoulder loop the more upright you will sit.

It is encouraged that the user keeps both arms on the outside of sling.

Start raising the hoist until the sling fits tightly around the back. Then tighten the adjustable belt with buckle around the waist. Ensure there is no discomfort.

Examination and Inspection of

General external wear

In normal use this is unavoidable and is shown by fluffiness to the surface fibres. This is harmless unless it becomes too extensive.

Local abrasion in the webbing straps

This may be caused by the passage of the webbing over sharp edges whilst under tension. Slight damage to the outer fibres may be considered safe but in severe cases this should lead to instant rejection of the product.

Cuts, holes or burns

These are dangerous and must lead to instant rejection of the product.

Examination of stitch patterns

Examine the stitch patterns for broken, worn, pulled or unravelled stitches. Reject any product where the thread has become subject to excessive abrasion or broken stitches.

Washing Instructions

Each sling can be laundered up to 80c using a mild soap solution and then rinse. Slings will air dry quickly in a warm room or in the normal manner. Do not place over a convector heater or steam pipes. Do not use Biological detergents. Do not tumble dry.

Safe Working Load

Parachute	200kg
Polyester	300kg
Polyester Mesh	300kg
Spacer Fabric	300kg